

# I.V. Stand (Intravenous pole) Instructions for Use

**Device Description:** An I.V. Stand is a stationary or mobile stand intended to hold infusion liquids, infusion accessories, and other medical devices.

**Intended Purpose:** The intended use of an I.V. Stand is to carry I.V. bags and pumps which helps assist ambulatory patients under I.V. therapy.

## Warnings, Cautions, and Contraindications:

### CAUTION:

- **MAXIMUM HOOK TOP LOAD: 8.8 lbs (4 kg)**
- **ADJUST UPPER POLE TO LOWEST POSSIBLE POSITION DURING TRANSIT.**
- **DO NOT OVERLOAD THIS STAND TO CREATE AN UNSTABLE SITUATION.**
- **NOT MRI COMPATIBLE UNLESS LABELED AS SUCH.**

### AMBULATORY SUPPORT WHEEL (OPTIONAL FEATURE)

- **MAXIMUM MASS OF USER- 300 lbs (136kg)**

### BRAKE CASTORS (OPTIONALFEATURE)

- **BRAKE CASTORS ENGAGED IN THE LOCKED POSITION INTRODUCES A RISK OF OVERBALANCE/TOPPLING UPON SUDDEN MOBILITY**

### TRANS-PAL® SUPPORT WHEEL: (OPTIONAL FEATURE) WHEN TRANS-PAL® IS ATTACHED TO A HEAD/FOOTBOARD THE FOLLOWING IS RECOMMENDED:

- **ADJUST BED HEIGHT TO LIFT IV STAND CASTORS OFF FLOOR BEFORE TRANSIT.**
- **KEEP HAND ON TRANS-PALSUPPORT WHEEL TO STABILIZE DURING TRANSIT.**
- **MOUNT PUMPS TO LOWEST POSSIBLE POSITION.**
- **USE CAUTION WHEN ROLLING OVER ELEVATOR THRESHOLDS AND OTHER UNEVEN SURFACES.**

**DO NOT DISMANTLE STAND FOR STORAGE OR TRANSPORT.  
Store indoors.**

## Assembly Instructions:

- **ASSEMBLY: FOLLOW PAGE ONE ILLUSTRATION NUMERAL SEQUENCE**

## TOOLS REQUIRED:

- **HEX KEY NO.1260 (SUPPLIED)**
- **WRENCH**

## Instructions for Use:

- **IV STAND HEIGHT ADJUSTMENT:**
  - **WHILEHOLDING INNER IV POLE, LOOSEN KNOB AND ADJUST TO DESIRED HEIGHT. RETIGHTEN KNOB.**
  - **GAS SPRING ADJUSTMENT: (OPTIONALFEATURE)**
    - **LOOSEN KNOB ON BASE POLE. LOWER OR RAISE SUPPORT WHEEL TO DESIRED HEIGHT. RETIGHTEN KNOB.**
  - **OXYGEN RINGS: (OPTIONAL FEATURE)**
    - **LOOSEN KNOB ON UPPER OXYGEN RING. LOWER OXYGEN CANISTER THROUGH RINGS TO REST ON CASTOR BASE. RETIGHTEN KNOB.**

## Connections with Other Devices:

Follow connecting instructions of the accessory manufacturer. Compatibility of accessories is not guaranteed.

## Cleaning:

It is recommended that the device is thoroughly wiped with a damp towel and a mild cleaning solution once a week or more often, as needed. Note: If IV solution spills on the base, clean immediately.

Cleaning solutions may be corrosive to this device. Follow hospital policies and procedures for the appropriate cleaning solution. Many hospitals prefer a mild solution of bleach, such as 10% bleach, a solution containing alcohol, or a solution of quaternary ammonium chloride. In diluted concentrations, these cleaning agents have no negative effects on Pryor Products IV stands, carts, or accessories.

For stainless steel, if rust occurs, use a Scotch-Brite™ pad to remove rust, and then use the above cleaning methods.


## Maintenance & Disposal:

Follow institutional policies and procedures, as well as local laws and regulations, when applicable.

# I.V. Stand (Intravenous pole) Instructions for Use

## Glossary:

Symbol	Name	Description
	Authorized representative in the European Community/ European Union	Indicates the authorized representative in the European Community/ European Union.
	Batch code	Indicates the manufacturer's batch code so that the batch or lot can be identified.
	Catalogue number	Indicates the manufacturer's catalogue number so that the medical device can be identified.
	Caution	To indicate that caution is necessary when operating the device or control close to where the symbol is placed, or to indicate that the current situation needs operator awareness or operator action in order to avoid undesirable consequences.
	CE Marking of Conformity	'CE marking of conformity' or 'CE marking' means a marking by which a manufacturer indicates that a device is in conformity with the applicable requirements set out in the EU Regulation 2017/745 and other applicable Union harmonisation legislation providing for its affixing;
	Consult instructions for use or consult electronic instructions for use	Indicates the need for the user to consult the instructions for use.
	Country of manufacture	To identify the country of manufacture of products. The date of manufacture is adjacent to this symbol.
	Manufacturer	Indicates the medical device manufacturer.
	Medical device	Indicates the item is a medical device.
	Packaging unit	To indicate the number of pieces in the package.
	Unique Device Identifier	Indicates a carrier that contains Unique Device Identifier information.

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