





# ELECTRICAL STRIP

## CAT. NO. 206

KIT PARTS LIST			
ITEM	DESCRIPTION	QTY	P/N
1	ELE. STRIP	1	MIS-1414
2	MOUNTING BAR	1	STL-1030
3	OXYGEN RINGS	1 SET	PLS-1403M
4	SQ. NUT	3	SCR-1256
5	SCREW	2	SCR-1609
6	CAP NUT	2	SCR-1484
7	SCREW	2	SCR-1302
8	PROTECTIVE COVER	1	MIS-1418M
9	SCREW	1	SCR-1667
10	KNOB	1	KNB-1259
11	TIE STRIP	1	MIS-1430
NOT SHOWN	HEX KEY	1	SCR-1260
	INSTRUCTION SHEET	1	LBL-3061



### ASSEMBLY INSTRUCTIONS:

• TOOLS REQUIRED: HEX KEY #1260 (SUPPLIED) PLIERS (NOT SUPPLIED)

#### • INSTALLATION OF OXYGEN RINGS TO A PRYOR I.V. STAND

- 1 REMOVE BOTH SCREWS THAT SECURE THE I.V. TOP POLE. REMOVE I.V. TOP POLE BY PULLING UP AND TWISTING LEFT TO RIGHT. (IF YOUR STAND HAS A SUPPORT WHEEL REMOVE SUPPORT WHEEL FIRST BY PULLING UP AND TWISTING.)
- 2 OXYGEN RINGS ARE BEING UTILIZED FOR ELECTRICAL STRIP MOUNTS. INSTALL OXYGEN RING ASSY (ITEM 3) AS SHOWN IN (FIG 1). NOTE: INSTALL OXYGEN RING WITH PROTECTIVE COVER AND KNOB LAST OVER THE BASE POLE. TIGHTEN SCREW (ITEM 5) VERY SECURELY AT THE DESIGNATED LOCATIONS LISTED IN (FIG 1) AGAINST THE BASE POLE.
- 3 INSTALL ELECTRICAL STRIP (ITEM 1) WITH CORD END DOWN, ALIGN THE MATING HOLE OF THE ELECTRICAL STRIP TO THE MATING SCREW HEADS (ITEM 7) LOCK INTO PLACE SECURELY.
- 4 SECURE ELECTRICAL STRIP (ITEM 1) BY WRAPPING STRIP TIE (ITEM 11) WITH SMOOTH SIDE FACING OUT AROUND THE MID SECTION OF ELECTRICAL STRIP (NOTE: TIE STRAP SHOULD BE LOCATED BETWEEN PLUG OUTLETS.) TIGHTEN SECURELY UNTIL THE STRAP TAKES THE SHAPE AND CONTOUR OF THE ELECTRICAL STRIP. (TRIM TIE AS REQUIRED). THE USE OF PLIERS (NOT SUPPLIED) MAY BE NEEDED FOR TIGHTENING.
- 6 REINSTALL SUPPORT WHEEL (IF YOUR STAND IS EQUIPPED) AND I.V. TOP POLE INTO CENTER HOLE. REINSERT AND TIGHTEN BOTH SCREWS SECURELY.

CLEANING: SOLUTIONS USED IN CONJUNCTION WITH THIS PRODUCT MAY BE CORROSIVE. IT IS RECOMMENDED TO WIPE DOWN IV STAND ONCE A WEEK MORE OFTEN IF REQUIRED WITH A DAMP TOWEL AND A MILD CLEANING SOLUTION. NOTE: IF IV SOLUTION SPILLS ON BASE, CLEAN IMMEDIATELY.